


The art of seeing
what others can't.



Handcrafted Loupes, Lights & Camera

admetec
We've got your back



We've got your back.





Admetec handcrafted, custom-designed loupes, lights, and Camera

We specialize in creating solutions that are exceptionally lightweight, convenient to use, and ergonomic - enabling you to focus on being effective while we focus on your comfort.

Our handcrafted loupes, lights, and Camera are custom designed for YOU. To do so, we take into consideration your unique facial structure, working distance, style, and personal preferences. User-centered product design combined with high-quality materials and components is the foundation of Admetec excellence.

Powerful lighting. Ultra-lightweight. Wireless.

Meet the Butterfly EVO, Admetec's cordless headlight series, featuring two models designed for optimal efficiency.

- 
Lightweight
 Ergonomically designed to rest lightly on your loupes
- 
Cordless
 Freedom of movement increases comfort and efficiency
- 
User-friendly controls
 Simple operation with tap controls for ON/OFF and light intensity
- 
Charging indicator display
 Always know your battery's charging status



Butterfly™ EVO

Butterfly combines effective lighting in an ultra-lightweight, easy-to-use, wireless configuration.

Intensity	Weight	Run time	Charge time
35k LUX	24 gr	2 hrs	1.25 hrs
Yellow filter	Batteries	Color temp	
External filter	3 pcs	5,750K	



Butterfly-S™ EVO

Butterfly-S provides stronger light intensity to enable work with higher magnifications.

Intensity	Weight	Run time	Charge time
55k LUX	29.5 gr	2 hrs	1.5 hrs
Yellow filter	Batteries	Color temp	
External filter	3 pcs	5,750K	



Convenient. Compact. Shockproof.

Your Admetec experience extends beyond working hours. What happens when you are not working with your loupes and light? Where do you put them? How do you protect and maintain them in good condition over time? If you work in more than one clinic, how do you take your loupes and light back and forth? That's why we designed convenient, compact, shockproof cases for each set of Admetec products.

Butterfly (all models)

3 batteries | Charging pad and cable | Loupes | User Manual



Orchid (all models)

2 PowerPacks | Charger | Loupes | User Manual

Loupes

Lanyard | Microfiber cloth | User Manual





We've got **your back.**

Admetec designs, develops, manufactures, and distributes handcrafted, custom-fitted magnifying loupes and LED headlamps as well as an integrated camera and headlight. We specialize in creating products that are exceptionally lightweight, convenient to use, and ergonomic – enabling you to focus on being effective while we focus on your comfort.

MM4000001

