

The art of seeing
what others can't.



Handcrafted Loupes, Lights & Camera

admetec
We've got your back



We've got your back.

Admetec handcrafted, custom-designed loupes, lights, and Camera

We specialize in creating solutions that are exceptionally lightweight, convenient to use, and ergonomic - enabling you to focus on being effective while we focus on your comfort.

Our handcrafted loupes, lights, and Camera are custom designed for YOU. To do so, we take into consideration your unique facial structure, working distance, style, and personal preferences. User-centered product design combined with high-quality materials and components is the foundation of Admetec excellence.

Galilean

Galilean loupes are extremely compact and lightweight while ensuring a wide field of view and excellent depth of field. Excellent for practitioners who want to see the big picture and don't necessarily need high magnification. Our Galilean loupes come in 3 magnifications.

Designed specifically for your working distance, to the centimeter.



2.5x

Working distance range
30-60 cm

Weight
51gr



2.7x

Working distance range
30-60 cm

Weight
43gr



3.2x

Working distance range
30-60 cm





Weight
43gr

What style makes your heart sing?

You're still YOU while you are at work. That's why we made sure to provide a wide variety of frames that allow you to express your unique personality and work in comfort.



Blues

-  Navy Blue
-  Bronze & Light Blue
-  Black
-  Pink







Indie

-  Champagne
-  Orange Lava
-  Midnight Blue
-  Black






Soul

-  Pink
-  Brown & Light Blue
-  Burgundy & Gold
-  Black







Jazz

-  Red & Black
-  Orange & Black
-  Blue & Black



Morriz of Sweden®

-  Red
-  Green
-  Blue
-  Black



Progear®

-  Red & Metallic Brown
-  Deep Green & White
-  White & Blue
-  Black



Bollé®

-  Gun Brown

Convenient. Compact. Shockproof.

Your Admetec experience extends beyond working hours. What happens when you are not working with your loupes and light? Where do you put them? How do you protect and maintain them in good condition over time? If you work in more than one clinic, how do you take your loupes and light back and forth? That's why we designed convenient, compact, shockproof cases for each set of Admetec products.

Butterfly (all models)

3 batteries | Charging pad and cable | Loupes | User Manual



Orchid (all models)

2 PowerPacks | Charger | Loupes | User Manual

Loupes

Lanyard | Microfiber cloth | User Manual





We've got **your back.**

Admetec designs, develops, manufactures, and distributes handcrafted, custom-fitted magnifying loupes and LED headlamps as well as an integrated camera and headlight. We specialize in creating products that are exceptionally lightweight, convenient to use, and ergonomic – enabling you to focus on being effective while we focus on your comfort.

MM4000001

